Guideline

Ligamys implants in the Magnetic Resonance (MR) Environment

Preservation in motion

*Building on our heritage
Moving technology forward
Step by step with our clinical partners
Towards a goal of preserving mobility*

Preservation in motion

As a Swiss company, Mathys is committed to this guiding principle and pursues a product portfolio with the goal of further developing traditional philosophies with respect to materials or design in order to address existing clinical challenges. This is reflected in our imagery: traditional Swiss activities in conjunction with continuously evolving sporting equipment.

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EN – Magnetic Resonance Imaging (MRI) safety information

The safety rating «MR Conditional» means that in a specific MR environment with defined application conditions a product demonstrably poses no hazards.

Notes on MRI safety

Non-clinical tests have shown Ligamys implants to be MR Conditional. A patient with such an implant can be safely examined in an MR system that satisfies the following conditions:

- Strength of the static magnetic field of 1.5 or 3.0 tesla, respectively
- Maximum spatial gradient field of 3 000 gauss/cm (30T/m)
- Maximum specific absorption rate specified for the MR system, whole-body average (SAR), of 4W/kg for 15 minutes of scanning

Note: The defined SAR must be strictly adhered to if the nearest part of the implant is closer than 30cm to the isocentre of the MR scanner. If the implant is in its entirety more than 30cm away from the isocentre of the MRI device, there is lower RF-induced heating of the tissue around the implants, and MRI scans in normal operating mode are permissible for all implants.

- Normal mode of the MR system. The effects of MRI procedures using MR systems and conditions above these levels have not been determined.
- The impact of local RF transmission coils has not been studied, and these are not recommended in the area of the implant.

Please note: Combination with other devices has not been evaluated for safety and compatibility in the MR environment. It has not been tested for RF heating in the MR environment. It is not known whether combination with other products is safe in the MR environment. Examination of a patient with such a product combination may result in patient injury. Therefore, combination of products has been rated as «MR unsafe».

Under the examination conditions defined above, it can be assumed that a Ligamys implant leads to a maximum temperature increase of less than 4.4°C after 15 minutes of continuous scanning. In non-clinical tests, the image artefact caused by the implant extends to about 55 mm from the Ligamys, if imaging is carried out using a gradient echo sequence and a 3.0-T MRI system. In non-clinical tests, the magnetically induced displacement force was tested. At 3.0T, a maximum deflection angle of less than 68° was found.

In non-clinical testing, the magnetically induced displacement torque was tested. No torque was measured.

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